PT

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EN

General instructions for use:
- Open the individual product package with all aseptic care.
- Verify all connections prior use.
- Purge the system. Make sure that there are no air-bubbles in the infusion system.
- To work by pressure, install the system on the infusion pump and connect to the patient.
- Place the drip chamber ±80 cm above the patient's insertion point.
- Open the flow regulator. The flow of the solution must be adjusted with the flow regulator at the beginning of the infusion and controlled during the procedure. To interrupt or stop the infusion, close the flow regulator.
- For rigid containers leave the vented spike air cap open.

Characteristics:
Sterile if pack un-opened and un-damaged. Non pyrogenic. Free of latex and phthalates. Standard connection to patient IV catheters. The maximum duration of use should be defined according to the recommendations of the Centre for Disease Control and Prevention (CDC), with replacement recommended every 24 hours. Prolonged use with various medications can also cause the components to weaken and consequently break and leak.

Intended use:
The infusion sets are single-use, sterile devices intended to channeling liquids for the purpose of infusion or administration into the body, working by gravity as well by pressure (< 2 bar) when connected to an infusion pump. A line that is suitable for gravity alone can, however, be indirectly connected to a pressure system. The intended users of these products are healthcare professionals (physicians and nurses) with knowledge and practical experience in handling these devices. It can be used for both adults and children.

Warnings:
Read Instructions for Use (IFU) before using the device. Do not use if the packaged is opened or damaged. Do not use if the label is illegible/ missing information. Do not use if protective caps are loose or missing. Verify all connections before use. If the product is damaged/out of specifications, do not use it and inform the manufacturer. The product is sterilized for single use only. Discard after single use. Please do not clean and do not reuse. Do not use solutions based on isopropyl alcohol. A reuse of the device may cause harmful infections, injury, or death. The infusion set is not suitable for infusion of blood and cell-containing blood products. The infusion of highly viscous solutions (ex: lipid emulsions and certain suspensions) can affect the flow rate. The product with the graduated chamber has no measurement function. Do not use after expiry date. If over it, please discard it. Store in a dry place and between 5 to 35°C. Discard all used devices according to the applicable legal regulations. The Summary of Product Characteristics can be consulted if in doubt about the compatibility of the medicine.

Performance:
Based on the instructions for use mentioned above and the capacity of the device to allow the ratio of each milliliter to be equivalent to 20 or 60 drops, the device presents a good performance.

Residual risks:
Although all risks are controlled and reduced to the minimum possible, there are always residual risks that can be limited, provided that these instructions and the practices of the clinical process are followed by healthcare professionals, namely: deterioration of the patient's state of health; infection, fever, embolism, underdosing, patient discomfort, toxic reactions.

Vigilance:
In case of serious incident related to the device, contact the manufacturer or the National Competent Authority.

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REF	Product Reference		Use-by date	(mm)	Tube Diameters	CE 2797	CE marking with Notified Body Number		Warning	15 µm	15 Micra pore size liquid filter	35°C 5°C	Temperature Limit	20 ml	Drops per milliliter	P Pressure	Compatible for use with pressure	2 STERILIZE	Do not resterilize	PHT DEHP SEM DEHP	Does not contain DEHP		Do not use if package is damaged		Keep way from sunlight	UDI	Unique Device Identifier		Manufacturer
	Date of manufacture	LOT	Batch code		Tube Dimensions	MD	Medical Device	eIFU indicator	Instructions for use in the QR code	TOTM	Plasticized in TOTM	AIR STOP	AirStop	STERILE EO	Sterilized using ethylene oxide	G Gravity	Compatible for use with gravity	VOL	Storage volume		Do not re-use	LATEX SEM LATEX	Does not contain natural rubber latex		Non-pyrogenic		Keep dry		Distributor

PMH - PRODUTOS MÉDICOS HOSPITALARES, SA
RUA GUINÉ-BISSAU LT.9, ZONA INDUSTRIAL DA MURTEIRA
2135-327 SAMORA CORREIA, PORTUGAL
Tel. +351 263 650 680 | WWW.PMH.PT

Implementatie datum /Implementation date /Data de implementação/Umsetzungsdatum/
Fecha de implementación/ Date de mise en oeuvre:28/10/2026