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
EN

Instructions for use:
- Remove aseptically, the Blood administration set from packaging.
- Verify all connections prior use.

Volume control administration set
-Close clamp and roller clamp. -Insert spike into blood bag. -Open the clamp. -Fill up calibrated chamber with necessary volume. -Close clamp. - Gently press and fill the blood chamber to ensure that the entire filter surface (inside and outside) is completely wetted. -Open clamp. -Purge the system. Make sure all air is removed by priming prior to use.

Set with blood chamber:
-Close roller clamp. -Remove the protective cap from spike. -Insert spike into blood bag. -Gently press and fill the blood chamber to ensure that the entire filter surface (inside and outside) is completely wetted. -Open the roller clamp. Make sure all air is removed by priming prior to use. -Close roller clamp.

Common Instructions to both systems
-Connect the system to patient. Open the roller clamp. The blood flow rate must be adjusted with a ramp, at the beginning and controlled during the transfusion. To interrupt or stop transfusion, close roller clamp.

Characteristics:
-Sterile if pack un-opened and un-damaged. Non pyrogenic. Free of latex and phthalates. -Standard connection to patient IV catheters. Limited exposure contact device (maximum contact duration of 24 hours).
It is recommended that no transfusion sets be used for more than 6 hours.
-The Centers for Disease Control, the Perfusion Nursing Society (USA) and other organizations publish useful guidelines for developing local guidelines. Consult local protocols.

Intended use:
The blood administration sets are single-use, sterile devices intended to channel blood, plasma or other components for the purpose of transfusion into the body.
The intended users of these products are healthcare professionals (physicians and nurses) with knowledge and practical experience in handling these devices. It can be used for both adults and children.

Warnings
- Read Instructions for Use (IFU) before using the device; -Do not use if the package is open or damaged.- do not use if the label is illegible/ missing information; -Do not use if protective caps are loose or missing. - Do not use if the filter is damaged or missing; -Verify all connections prior use; If the product is damaged/out of specifications, do not use it and inform the manufacturer; -The product is sterilized for single use only. Discard after single use. -Please do not clean and do not reuse. Do not use solutions based on isopropyl alcohol. A reuse of the device may cause harmful infections, injury or death. **-The product with the graduated chamber has no measurement function**. -Do not use after expiry date. If over it, please discard it; - Store in a dry place and between 5 to 35°C; - Cleaning solutions and disinfectants may deteriorate the materials used for bloodlines and should not be used.
Used products are subject to disposal and destruction procedures in accordance with current legislation in the country in which the Bloodlines are used. If necessary, waste will be disposed of by an accredited organization in accordance with applicable legislation.

Performance:
Based on the instructions for use mentioned above and the capacity of the device to allow the ratio of each milliliter to be equivalent to 20 drops, the device presents a good performance.

Residual Risks:
Although all risks are under control and as low as possible, there are always residual risks that can be limited as long as these instructions and clinical/transfusion process practices are followed by healthcare professionals, namely: malfunction due to leakage or occlusion, infection, fever, embolism, haemolysis, patient discomfort, toxic reactions and deterioration of the patient's health (anaemia, inflammatory reactions).

Vigilance
In case of serious incident related to the device, contact the manufacturer or the National Competent Authority.
If additional operating instructions are required, these can be requested from the manufacturer or can be obtained from the PMH homepage: https://www.pmh.pt/familia200info/

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REF	Product Reference		Use-by date	(mm)	Tube Diameters	CE 2797	CE marking with Notified Body Number		Warning	15 µm	15 Micra pore size liquid filter	35°C 5°C	Temperature Limit	20 ml	Drops per milliliter	P Pressure P	Compatible for use with pressure	2 STERILIZE	Do not resterilize	PHT DEHP SEM DEHP	Does not contain DEHP		Do not use if package is damaged		Keep way from sunlight	UDI	Unique Device Identifier		Manufacturer
	Date of manufacture	LOT	Batch code		Tube Dimensions	MD	Medical Device	eIFU indicator	Instructions for use in the QR code	TOTM	Plasticized in TOTM	AIR STOP	AirStop	STERILE EO	Sterilized using ethylene oxide	G Gravity G	Compatible for use with gravity	VOL	Storage volume	2	Do not re-use	LATEX SEM LATEX	Does not contain natural rubber latex		Non-pyrogenic		Keep dry		Distributor

PMH - PRODUTOS MÉDICOS HOSPITALARES, SA
RUA GUINÉ-BISSAU LT.9, ZONA INDUSTRIAL DA MURTEIRA
2135-327 SAMORA CORREIA, PORTUGAL
Tel. +351 263 650 680 | WWW.PMH.PT

Implementatie datum /Implementation date /Data de implementação/Umsetzungsdatum/
Fecha de implementación/ Date de mise en oeuvre: 28/10/2026