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EN

Instructions for use:
- Open the individual product package with all aseptic care.
- Verify all connections prior use.
- Purge the system. Make sure that all air bubbles are removed by priming prior to use.
- Connect to the necessary infusion set and/or to the catheter.
- This extension set is intended to work by gravity or connected to a pressure pump.
- Discard all used devices according to the applicable legal regulation. Specially marked containers should be used if necessary.

Characteristics:
-Sterile if pack unopened and undamaged
-Free of latex and DEHP
-Non-pyrogenic
-Standard connectors.
-Medical device for single use.
- The maximum duration of use should be defined according to the recommendations of the Centre for Disease Control and Prevention (CDC), with replacement recommended every 24 hours. Prolonged use with various medications can also weaken the components, causing them to break and leak.

Intended use:
The extension sets are single-use, sterile devices intended to channel liquids for the purpose of infusion or administration into the body, working by gravity as well by pressure when connected to an infusion pump, according to the labelling of the device. A line that is suitable for gravity alone can, however, be indirectly connected to a pressure system. For contrast administration, only use devices labelled accordingly or with a compatible maximum pressure indication. Some lines can be used for liquids withdraw.
The intended users of these products are healthcare professionals (physicians and nurses) with knowledge and practical experience in handling these devices. It can be used for both adults and children.

Warnings:
Read Instructions for Use (IFU) before using the device. Do not use if the packaging is open or damaged. Do not use if the label is illegible/ missing information. Do not use if protective caps are loose or missing. Verify all connections before use. Do not apply excessive force when connecting or disconnecting the device, as this may damage or fracture the device and result in leakage. Use only standard, compatible connectors in accordance with ISO 80369-7. If the product is damaged/out of specifications, do not use it and inform the manufacturer. The device is sterilized for single use only. Discard after single use. Do not clean or reuse. This product is not intended for reuse. A reuse of the device may cause harmful infections, injury, or death. Do not use after expiry date. If over it, please discard it. This device is not suitable for isopropyl alcohol-based solutions. Store in a dry place and between 5 to 35°C. Discard all used devices according to the applicable legal regulations. The Summary of Product Characteristics can be consulted if in doubt about the compatibility of the medicine.

Performance:
Based on the instructions for use mentioned above and the capacity of the device to allow visualization of a continuous flow without leaks or obstruction, the device presents a good performance.

Residual risks:
Although all risks are controlled and reduced to the minimum possible, there are always residual risks that can be limited, provided that these instructions and the practices of the clinical process are followed by healthcare professionals, namely: deterioration of the patient's state of health; infection, fever, embolism, underdosing, patient discomfort, toxic reactions.

Vigilance:
In case of serious incident related to the device, contact the manufacturer or the National Competent Authority.
If further Instructions for Use are required, these can be requested from the manufacturer or can be obtained from the PMH homepage: https://www.pmh.pt/familia300info/

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REF	Product Reference		Use-by date		Tube Diameters	CE 2797	CE marking with Notified Body Number		Warning	15 µm	15 Micra pore size liquid filter	35°C 5°C	Temperature Limit	20 ml	Drops per milliliter	P Pressure	Compatible for use with pressure	2 STERILIZE	Do not resterilize	PHT DEHP SEM DEHP	Does not contain DEHP		Do not use if package is damaged		Keep way from sunlight	UDI	Unique Device Identifier		Manufacturer
	Date of manufacture	LOT	Batch code		Tube Dimensions	MD	Medical Device	IFU indicator	Instructions for use in the QR code	TOTM	Plasticized in TOTM	AIR STOP	AirStop	STERILE EO	Sterile barrier. Sterilized using ethylene oxide	G Gravity	Compatible for use with gravity	VOL	Storage volume		Do not re-use	LATEX SEM LATEX	Does not contain natural rubber latex		Non-pyrogenic		Keep dry		Distributor

PMH - PRODUTOS MÉDICOS HOSPITALARES, SA
RUA GUINÉ-BISSAU LT.9, ZONA INDUSTRIAL DA MURTEIRA
2135-327 SAMORA CORREIA, PORTUGAL
Tel. +351 263 650 680 | WWW.PMH.PT

Implementatie datum /Implementation date /Data de implementação/Umsetzungsdatum/
Fecha de implementación/ Date de mise en oeuvre: 28/10/2026