


**DESCRIPTION**

- Double channel spike with air inlet with filter
- Tube 15 + 5 cm
- Clamps
- Graduated chamber / air inlet with filter / cap with injection site
- Connector

ALTERNATIVE	REFERENCE	UNITS/BOX
Chamber 120 ml	109.0745 / SMR120	50
Chamber 150 ml	109.0746 / SMR150	50

STERILE EO

Ref.: 109.0925 Units/Box: 50



**DESCRIPTION**

- Simultaneous administration solutions set
- 2 Spikes
- Tube 15 + 160 cm
- Clamps
- Precision roller clamp
- Graduated chamber 120 ml / 20 dps/ml / air inlet with filter / filter 32 cm<sup>2</sup> / valve
- Injection site latex free
- Luer lock

STERILE EO




**DESCRIPTION**

- Double channel spike with air inlet with filter
- Tube opaque 15 + 160 cm
- Precision roller clamps
- Graduated chamber microdrip 120 ml / 60dps/ml / opaque / air inlet with filter / cap with injection site
- Luer lock

ALTERNATIVE	REFERENCE	UNITS/BOX
Tube opaque black	109.0227 / SM120LOPA	50
Tube opaque orange	109.0853 / SM120LOPALAR	50


STERILE EO



**DESCRIPTION**

- Double channel spike with air inlet with filter
- Tube 15 + 60 + 100 cm
- Precision roller clamps
- Graduated chamber microdrip / 60 dps/ml / air inlet with filter / cap with injection site
- Flow control device – STAT2
- Injection site
- Luer lock


ALTERNATIVE	REFERENCE	UNITS/BOX
Chamber 120 ml	109.0242 / STAT2SM120B	50
Chamber 150 ml	109.0780 / STAT2SM150B	50



**DESCRIPTION**

- Double channel spike with air inlet with filter
- Tube 15 + 160 cm
- Precision roller clamps
- Graduated chamber microdrip 120 ml / 60 dps/ml / air inlet with filter / cap with injection site
- Injection site
- Luer lock

ALTERNATIVE	REFERENCE	UNITS/BOX
Chamber microdrip with filter 15µm	109.0864 / SM120BF	50
Chamber microdrip without filter	109.0223 / SM120B	50




**DESCRIPTION**

- Double channel spike with air inlet with filter
- Tube 15 + 160 cm
- Precision roller clamps
- Graduated chamber microdrip 150 ml / 60 dps/ml / air inlet with filter / cap with injection site
- Injection site
- Luer lock

ALTERNATIVE	REFERENCE	UNITS/BOX
Chamber microdrip with filter 15µm	109.0232 / SM150BF	50
Chamber microdrip without filter	109.0231 / SM150B	50






**DESCRIPTION**

- Double channel spike with air inlet with filter
- Tube 15 + 140 + 20 cm
- Precision roller clamps
- Graduated chamber microdrip 60 dps/ml / filter 15µm / air inlet with filter / cap with injection site
- Y injection site
- Luer lock

ALTERNATIVE	REFERENCE	UNITS/BOX
Chamber microdrip 120 ml	109.0229 / SM120YF	50
Chamber microdrip 150 ml	109.0236 / SM150YF	50

STERILE EO  



**DESCRIPTION**


- Double channel spike with air inlet with filter
- Tube 15 + 160 cm
- Precision roller clamps
- Graduated chamber microdrip 120 ml / 60 dps/ml / air inlet with filter / cap with injection site latex free / valve
- Injection site latex free
- Luer lock

ALTERNATIVE	REFERENCE	UNITS/BOX
Chamber microdrip with filter 15µm	109.0224 / SM120BFV	50
Chamber microdrip without filter	109.0225 / SM120BV	50

STERILE EO 

Ref.: 109.0228

Units/Box: 50




**DESCRIPTION**

**SM120VYF**

- Double channel spike with air inlet with filter
- Tube 15 + 140 + 20 cm
- Precision roller clamps
- Graduated chamber microdrip 120 ml / 60 dps/ml / air inlet with filter / cap with injection site latex free / valve / filter 15µm
- Y injection site latex free
- Luer lock

STERILE EO 



**DESCRIPTION**

**SM150BFV**

- Double channel spike with air inlet with filter
- Tube 15 + 160 cm
- Precision roller clamps
- Graduated chamber microdrip 150 ml / 60 dps/ml / air inlet with filter / cap with injection site latex free / filter 15µm / valve
- Injection site latex free
- Luer lock

STERILE EO




**DESCRIPTION**

**SM150BFVO**

- Double channel spike with air inlet with filter
- Tube opaque black 15 + 160 cm
- Precision roller clamps
- Graduated chamber microdrip 150 ml / 60 dps/ml / air inlet with filter / cap with injection site latex free / opaque / filter 15µm / valve
- Injection site latex free
- Luer lock

STERILE EO




**DESCRIPTION**

- Double channel spike with air inlet with filter
- Tube 15 + 140 + 20 cm
- Precision roller clamps
- Graduated chamber microdrip 150 ml / 60 dps/ml / air inlet with filter / cap with injection site latex free / valve / filter 15µm
- Y injection site latex free with microClave® - valve for needle free injection/aspiration (with disinfection)
- Luer lock

STERILE EO



**DESCRIPTION**

- Double channel spike with air inlet with filter
- Tube 15 + 140 + 20 cm
- Precision roller clamps
- Graduated chamber microdrip 150 ml / 60 dps/ml / air inlet with filter / cap with injection site latex free / valve / filter 15µm
- Y injection site latex free with valve for needle free injection/aspiration (with disinfection)
- Luer lock

STERILE EO

