

**DESCRIPTION****QUIMIO2PUMP**

- 2 Way cytosstatics administration set for pump
- Double channel spike with air inlet with filter
- Tube 12 + 10 cm
- Clamps
- 3 Way connection with 2 VAR1 - anti-reflux valve for needle free injection

**DESCRIPTION****QUIMIO2PUMPOPA**

- 2 Way cytosstatics administration set for pump
- Double channel spike with air inlet with filter
- Tube opaque orange 12 + 10 cm
- Clamps
- 3 Way connection with 2 VAR1 - anti-reflux valve for needle free injection

**DESCRIPTION****QUIMIO4PUMP**

- 4 Way cytosstatics administration set for pump
- Double channel spike with air inlet with filter
- Tube 12 + 5 + 10 cm
- Clamps
- 2 Three way connections with 4 VAR1 - anti-reflux valve for needle free injection

**DESCRIPTION****QUIMIO4PUMPOPA**

- 4 Way cytosstatics administration set for pump
- Double channel spike with air inlet with filter
- Tube opaque orange 12 + 5 + 10 cm
- Clamps
- 2 Three way connections with 4 VAR1 - anti-reflux valve for needle free injection





DESCRIPTION

- Cytostatics administration set for pump
- Double channel spike with air inlet with filter
- Tube 30 + 10 cm
- Clamps
- Y injection site with microClave® - valve for needle free injection/aspiration (with disinfection)


STERILE EO



DESCRIPTION

- Cytostatics administration set for pump
- Double channel spike with air inlet with filter
- Tube 30 + 10 cm
- Clamps
- Y injection site with valve for needle free injection /aspiration (with disinfection)

STERILE EO



DESCRIPTION

- 2 Way cytotstatics administration set for pump
- Double channel spike with air inlet with filter
- Tube 15 + 10 cm
- Clamps
- 3 Way connection with 2 valve for needle free injection/aspiration (with disinfection)

STERILE EO



DESCRIPTION

- 4 Way cytotstatics administration set for pump
- Double channel spike with air inlet with filter
- Tube 15 + 2,5 + 10 cm
- Clamps
- 2 Three way connections with 4 valve for needle free injection/aspiration (with disinfection)

STERILE EO

**DESCRIPTION****QUIMIO2**

- 2 Way cytosstatics administration set
- Double channel spike with air inlet with filter
- Tube 12 + 3 + 160 cm
- Clamp
- 3 Way connection with 2 VAR 1 - anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Precision roller clamp
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



NO DEHP

**DESCRIPTION****QUIMIO2OPA**

- 2 Way cytosstatics administration set
- Double channel spike with air inlet with filter
- Tube opaque orange 12 + 3 + 160 cm
- Clamp
- 3 Way connection with 2 VAR1 – anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / opaque / filter 15µm
- Precision roller clamp
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



NO DEHP

**DESCRIPTION****QUIMIO2VAR2**

- 2 Way cytosstatics administration set
- Double channel spike with air inlet with filter
- Tube 12 + 3 + 160 + 50 cm
- Precision roller clamps
- 3 Way connection with 2 VAR1 - anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Y injection site with VAR2 - anti-reflux valve for needle free injection/aspiration
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



NO DEHP

**DESCRIPTION****QUIMIO2TR3**

- 2 Way cytosstatics administration set
- Double channel spike with air inlet with filter
- Tube 12 + 5 + 160 + 50 cm
- Precision roller clamps
- 3 Way connection with 2 VAR1– anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Luer lock / 3 way stopcock
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



NO DEHP



DESCRIPTION

QUIMIO4

- 4 Way cystostatics administration set
- Double channel spike with air inlet with filter
- Tube 12 + 5 + 3 + 160 cm
- Clamp
- 2 Three way connections with 4 VAR1- anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Precision roller clamp
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



DESCRIPTION

QUIMIO4OPA

- 4 Way cystostatics administration set
- Double channel spike with air inlet with filter
- Tube opaque orange 12 + 5 + 3 + 160 cm
- Clamp
- 2 Three way connections with 4 VAR1 - anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / opaque / filter 15µm
- Precision roller clamp
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



DESCRIPTION

QUIMIO4TR3

- 4 Way cystostatics administration set
- Double channel spike with air inlet with filter
- Tube 12 + 5 + 3 + 160 + 50 cm
- Precision roller clamps
- 2 Three way connections with 4 VAR1- anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Luer lock / 3 way stopcock
- Luer lock mobile
- Cap with automatic purge filter


STERILE EO



DESCRIPTION

- Cystostatics administration set
- Double channel spike PVC free with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Tube 145 + 15 cm
- Precision roller clamp
- Hydrophilic filter 0,2µm
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO


**DESCRIPTION**

- Cytostatics administration set
- Double channel spike PVC free with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Tube 200 + 15 + 14 cm
- Precision roller clamp
- Hydrophilic filter 0,2µm
- Luer lock / 3 way stopcock
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



NO DEHP


**DESCRIPTION**

- Cytostatics administration set
- Double channel spike PVC free with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Tube 145 + 15 + 20 cm
- Precision roller clamp
- Hydrophilic filter 0,2µm
- Y injection site
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



NO DEHP


**DESCRIPTION**

- Double channel spike PVC free with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking
- Tube 140 + 80 + 20 cm
- Flow control device - EASYDROP
- Clamp
- Y injection site with microClave® - anti-reflux valve for needle free injection/aspiration (with disinfection)
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO



NO DEHP



DESCRIPTION

- Pump administration set
- Double channel spike with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking
- Special tube for pump 150 cm
- Precision roller clamp
- Luer lock

ALTERNATIVE	REFERENCE	UNITS/BOX
Chamber with filter 15µm	110.0171 / SS200FLBOMBA	200
Chamber without filter 15µm	110.0216 / SS200LBOMBA	200

STERILE EO NO DEHP

Ref.: 110.0241

Units./Box: 50



DESCRIPTION

SS120LBOMBA

- Pump administration set
- Double channel spike with air inlet with filter
- Tube 15 cm
- Precision roller clamp
- Graduated Chamber 120 ml / 20 dps/ml / air inlet with filter / cap with injection site
- Special tube for pump 160 cm
- Precision roller clamp orange
- Luer lock

STERILE EO NO DEHP

Ref.: 110.0226

Units./Box: 50



DESCRIPTION


SM120LBOMBA

- Pump administration set
- Double channel spike with air inlet with filter
- Tube 15 cm
- Precision roller clamp
- Graduated Chamber 120 ml / 60 dps/ml / air inlet with filter / cap with injection site
- Special tube for pump 160 cm
- Precision roller clamp orange
- Luer lock

STERILE EO NO DEHP

Ref.: 110.1688

Units./Box: 50



DESCRIPTION

- Double channel spike with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Special tube for pump 160 + 30 cm
- Precision roller clamp orange
- Y injection site with microClave® - anti-reflux valve for needle free injection/aspiration (with disinfection)
- Luer lock
- Cap with automatic purge filter

STERILE EO NO DEHP