**DESCRIPTION**

- 2 Way cytosstatics administration set
- Double channel spike with air inlet with filter
- Tube POE* 12 + 5 + 160 + 50 cm
- Precision roller clamps
- 3 Way connection with 2 VAR1 - anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Luer lock / 3 way stopcock
- Luer lock
- Cap with automatic purge filter

*POE - Polyolefin

**DESCRIPTION**

- 2 Way cytosstatics administration set
- Double channel spike with air inlet with filter
- Tube opaque orange POE* 12 + 5 + 160 + 50 cm
- Precision roller clamps
- 3 Way connection with 2 VAR1 - anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / opaque / filter 15µm
- Luer lock / 3 way stopcock
- Luer lock
- Cap with automatic purge filter

*POE - Polyolefin

**DESCRIPTION**

- 4 Way cytosstatics administration set
- Double channel spike with air inlet with filter
- Tube POE* 12 + 3 + 5 + 160 + 50 cm
- Precision roller clamp
- 2 Three way connections with 4 VAR1 - anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Luer lock / 3 way stopcock
- Luer lock
- Cap with automatic purge filter

*POE - Polyolefin






DESCRIPTION

- 2 Way cytotostatics administration set
- Double channel spike with air inlet with filter
- Tube opaque orange POE* 12 + 3 + 5 + 160 + 50 cm
- Precision roller clamps
- 2 Three way connections with 4 VAR1 - anti-reflux valve for needle free injection
- Chamber PVC free 20 dps/ml / solution marking / opaque / filter 15µm
- Luer lock / 3 way stopcock
- Luer lock
- Cap with automatic purge filter

*POE - Polyolefin



DESCRIPTION

- Cystostatics administration set
- Double channel spike with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Tube POE* 150 cm
- Precision roller clamp
- Luer lock

*POE - Polyolefin

