


DESCRIPTION

- Double channel spike with air inlet with filter
- Tube 15 + 140 + 20 cm
- Precision roller clamps
- Graduated chamber microdrip 150 ml / 60 dps/ml / air inlet with filter / cap with injection site latex free / valve / filter 15µm
- Y injection site latex free with valve for needle free injection/aspiration (with disinfection)
- Luer lock

STERILE	EO
---------	----






DESCRIPTION

PSIUR2P

- High flow irrigation set
- 2 Double channel spikes with air inlet with filter
- Tube 60 + 150 cm
- 3 Clamps
- Y connection
- Connector

STERILE EO




DESCRIPTION

PSART2P

- High flow irrigation set
- 2 Spikes
- Tube 60 + 150 cm
- 3 Clamps
- Y connection
- Luer lock mobile

STERILE EO




DESCRIPTION

PSLV2P

- High flow irrigation set
- 2 Spikes
- Tube 20 + 170 cm
- 2 Clamps
- Y connection
- Roller clamp
- Connector with tube SEBS* 10 cm

* SEBS - Copolymers Hydrogenated Polybutadiene

STERILE EO




DESCRIPTION

PSLV1PC

- High flow irrigation set
- Spike
- Tube 40 + 150 cm
- Clamp
- Chamber PVC
- Roller clamp
- Luer lock mobile


STERILE EO



**DESCRIPTION****PSART2PC**

- High flow irrigation set
- 2 Spikes
- Tube 20 + 150 cm
- 3 Clamps
- Y connection
- Chamber PVC
- Luer lock mobile


STERILE	EO	2	LATEX
---------	----	---	-------

**DESCRIPTION****PSLV2PC**

- High flow irrigation set
- 2 Spikes
- Tube 20 + 170 cm
- 2 Clamps
- Y connection
- Chamber PVC
- Roller clamp
- Connector with tube SEBS* 10 cm

* SEBS - Copolymers Hydrogenated Polybutadiene

STERILE	EO	2	LATEX
---------	----	---	-------

**DESCRIPTION****PSPO**


- High flow irrigation set
- Spike
- Chamber PVC
- Tube 320 cm
- Roller clamp
- Connector

STERILE	EO	2	LATEX
---------	----	---	-------

**DESCRIPTION****2PSPO**

- High flow irrigation set
- 2 Spikes
- 2 Chambers PVC
- Tube 60 + 260 cm
- 3 Roller clamps
- Y Connection
- Connector

STERILE	EO	2	LATEX
---------	----	---	-------




DESCRIPTION

SS200FHLF

- Double channel spike with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking / hydrophilic filter 5µm
- Tube 150 cm
- Precision roller clamp
- Luer lock
- Cap with automatic purge filter

STERILE EO




DESCRIPTION

- Double channel spike with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking / hydrophilic filter 5µm
- Tube 140 + 20 cm
- Precision roller clamp
- Y injection site
- Luer lock
- Cap with automatic purge filter

ALTERNATIVE	REFERENCE	UNITS/BOX
Luer lock	102.0052 / SS200FHYF	200
Luer lock mobile	102.0058 / SS200FHYLF	200

STERILE EO



DESCRIPTION

SS200FHYVAR3

- Double channel spike with air inlet with filter
- Chamber PVC free 20 dps/ml / solution marking / hydrophilic filter 5µm
- Tube 140 + 20 cm
- Precision roller clamp
- Y injection site with VAR 3 - valve for needle free injection / aspiration (with disinfection)
- Luer lock
- Cap with automatic purge filter

STERILE EO

