



DESCRIPTION

IV²QUIMIO2

- 2 Way cytotostatics administration set
- Double channel spike with air inlet with filter
- Tube 12 + 4 + 3 + 160 cm
- Clamp
- 3 Way connection with 2 VAR1- anti-reflux valve for needle free injection
- **IV² flow** – equal pressure chamber
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Precision roller clamp
- Luer lock mobile
- Cap with automatic purge filter



DESCRIPTION

IV²QUIMIO2OPA

- 2 Way cytotostatics administration set
- Double channel spike with air inlet with filter
- Tube opaque orange 12 + 4 + 3 + 160 cm
- Clamp
- 3 Way connection with 2 VAR1 - anti-reflux valve for needle free injection
- **IV² flow** - equal pressure chamber
- Chamber PVC free 20 dps/ml / solution marking / opaque / filter 15µm
- Precision roller clamp
- Luer lock mobile
- Cap with automatic purge filter



DESCRIPTION

IV²QUIMIO4

- 4 Way cytotostatics administration set
- Double channel spike with air inlet with filter
- Tube 12 + 5 + 4 + 3 + 160 cm
- Clamp
- 2 Three way connections with 4 VAR1 - anti-reflux valve for needle free injection
- **IV² flow** - equal pressure chamber
- Chamber PVC free 20 dps/ml / solution marking / filter 15µm
- Precision roller clamp
- Luer lock mobile
- Cap with automatic purge filter





DESCRIPTION

IV²QUIMIO40PA

- 4 Way cystostatics administration set
- Double channel spike with air inlet with filter
- Tube opaque orange 12 + 5 + 4 + 3 + 160 cm
- Clamp
- 2 Three way connections with 4 VAR1 - anti-reflux valve for needle free injection
- **IV² flow** - equal pressure chamber
- Chamber PVC free 20 dps/ml / solution marking / opaque / filter 15µm
- Precision roller clamp
- Luer lock mobile
- Cap with automatic purge filter

STERILE EO   NO DEHP