
**DESCRIPTION****PR**


- Ruler for central venous pressure monitoring set

**DESCRIPTION****PPVC**

- Central venous pressure monitoring set without ruler
- Luer slip / cap with automatic purge filter
- Tube 150 + 60 cm
- Luer lock / 3 way stopcock
- Luer lock


**DESCRIPTION****PPVC R**

- System with ruler for central venous pressure monitoring set
- Ruler
- Luer slip / cap with purge filter
- Tube 150 + 60 cm
- Luer lock / 3 way stopcock
- Luer lock

**DESCRIPTION****PPVCRSS**

- System with ruler for central venous pressure monitoring set
- Ruler
- Luer slip / cap with purge filter
- Tube 150 + 60 cm
- Luer lock / 3 way stopcock
- Luer lock
- IV Sets (108.0201 / SS200FL)






**DESCRIPTION**

**PPIC**

- System with ruler for cerebral spinal fluid pressure monitoring set
- Tube 100 cm
- Luer lock mobile

STERILE EO



**DESCRIPTION**

**PKMA01**

- Arterial monitoring set
- Chamber PVC with spike
- Tube 160 cm
- Precision roller clamp
- Luer lock / TRANSDUCER
- Luer lock female
- Tube with red line 150 cm
- Luer lock mobile / 3 way stopcock red / cap with injection site
- Luer lock female
- Tube with red line 30 cm
- Luer lock mobile

STERILE EO



**DESCRIPTION**

**PKMA01MC**

- Arterial monitoring set
- Chamber PVC with spike
- Tube 160 cm
- Precision roller clamp
- Luer lock / TRANSDUCER
- Luer lock female
- Tube with red line 150 cm
- Luer lock mobile / 3 way stopcock red / microClave® - anti-reflux valve for needle free injection/aspiration (with disinfection)
- Luer lock female
- Tube with red line 30 cm
- Luer lock mobile


STERILE EO



**DESCRIPTION****PKMV01**

- Venous monitoring set
- Chamber PVC with spike
- Tube 160 cm
- Precision roller clamp
- Luer lock / TRANSDUCER
- Luer lock female
- Tube with blue line 150 cm
- Luer lock mobile / 3 way stopcock blue / cap with injection site
- Luer lock female
- Tube with blue line 30 cm
- Luer lock mobile

STERILE EO


**DESCRIPTION****PKMA02**

- Arterial/venous monitoring set
- Chamber PVC with spike
- Tube 160 cm
- Precision roller clamp
- Y connection:

**a)** tube 30 cm / clamp / luer lock / TRANSDUCER / luer lock female / tube with blue line 150 cm / luer lock mobile / 3 way stopcock blue / cap with injection site / luer lock female / tube with blue line 30 cm / luer lock mobile

**b)** tube 30 cm / clamp / luer lock / TRANSDUCER / luer lock female / tube with red line 150 cm / luer lock mobile / 2 three way stopcocks red / cap with injection site / luer lock female / tube with red line 30 cm / luer lock mobile

STERILE EO

**DESCRIPTION****PKMA02MC**

- Arterial/venous monitoring set
- Chamber PVC with spike
- Tube 160 cm
- Precision roller clamp
- Y connection:

**a)** tube 30 cm / clamp / luer lock / TRANSDUCER / luer lock female / tube with blue line 150 cm / luer lock mobile / 3 way stopcock blue / microClave® / luer lock female / tube with blue line 30 cm / luer lock mobile

**b)** tube 30 cm / clamp / luer lock / TRANSDUCER / luer lock female / tube with red line 150 cm / luer lock mobile / 2 three way stopcocks red / microClave® / luer lock female / tube with red line 30 cm / luer lock mobile microClave® - anti-reflux valve for needle free injection/aspiration (with disinfection)

STERILE EO

